DESCRIPTION: 30 mL Sensitive Eyes Daily Cleaner F/E Carton Canada

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BAUSCH + LOMB
Sensitive Eyes®
Daily Cleaner

BAUSCH + LOMB
Sensitive Eyes®
Daily Cleaner

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DAILY CLEANER

DIRECTIONS:

- After removing the lens, place it in the palm of your hand, and place three drops of Sensitive Eyes Daily Cleaner on each lens surface.
- Rub the lens for 20 seconds. Be sure to clean each side of the lens.
- Rinse the lens thoroughly with Bausch & Lomb Sensitive Eyes® Saline Plus Solution.
- Disinfect using the procedure which your eye care practitioner has recommended.

PRECAUTIONS:

- Keep this bottle tightly closed when not in use.
- Store at room temperature.
- Use before the expiration date marked on the carton and bottle.

SEE INSIDE OF CARTON FOR COMPLETE INSTRUCTIONS AND IMPORTANT SAFETY INFORMATION.

CONTENTS: A sterile, isotonic, buffered solution that contains hydroxypropyl methylcellulose, polyamine, sodium borate and sodium chloride, preserved with edetate disodium (0.5%) and sorbic acid (0.2%).

If you are allergic to any ingredient in this product, do not use. NOT INTENDED FOR USE DIRECTLY IN THE EYE.

621102

Removes film, deposits and debris

For use with all soft contact lens disinfection systems

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CUT HERE

ACTIONS: Dissolve and eliminate the impurities, deposits and debris from your lenses.

INDICATIONS: Bausch & Lomb Sensitive Eyes Daily Cleaner is indicated for use in cleaning soft contact lenses, in conjunction with thermal disinfection systems. This cleaner may be used with standard wear soft contact lenses as often as daily or as recommended by your eye care practitioner.

CONTRAINDICATIONS: If you are allergic to any ingredient in Bausch & Lomb Sensitive Eyes Daily Cleaner, do not use this product.

PRECAUTIONS:

- Use before the expiration date marked on the bottle and carton.
- Use before the lot number is printed on the bottle and carton.
- Replace cap after each use.
- Do not use if the seal is broken.
- Not intended for use directly in the eye.

ADVERSE REACTIONS: The following problems may occur eyes sting, burn, itch, or become watery. If you experience any of these symptoms, discontinue use immediately and consult your eye care practitioner.

IMMEDIATELY REMOVE YOUR LENSES.

- If the discomfort or problem stops, then look closely at the lens.
- If the discomfort or problem persists, then look closely at the lens.
- If the lens has dirt, an eyelash, or other foreign material on it, remove the lens by the frame and clean it thoroughly with clean saline, and disinfect the lens, then reinsert it.
- If the problem continues, IMMEDIATELY remove the lens and consult your eye care practitioner.

DIRECTIONS:

- After removing the lens, place it in the palm of your hand, and place three drops of Sensitive Eyes Daily Cleaner on each lens surface.
- Rub the lens for 20 seconds while in the palm of your hand using your thumb and index finger.
- Rinse the lens thoroughly by holding the lens in a steady stream of water for 20 seconds.
- Disinfect using the procedure which your practitioner has recommended.

NOTE: This product is not intended to remove protein deposits.

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©/TM are trademarks of Bausch & Lomb Incorporated or its affiliates. MD/UMC sont des marques de commerce de Bausch & Lomb Incorporated ou de ses filiales.
©/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.